


# Periodontal B

## PROPHY SET UP 15 INS (Recommended)

1. 1018 3EZ EXPLORERS OUD11/12
2. 1132 3EZ COLOR PROBE EXP23/CP12
3. 1201 6EZ SCALERS/CURETTES H6/7
4. 1208 6EZ SCALERS/CURETTES 204S
5. 1416 6EZ SCALERS/CURETTES COL13/14
6. 1601 6EZ STANDARD GRACEY GR1/2
7. 1604 6EZ STANDARD GRACEY GR7/8
8. 1606 6EZ STANDARD GRACEY GR11/12
9. 1607 6EZ STANDARD GRACEY GR13/14
10. 6517 FORCEP LONDON-COLLEGE
11. 6911 6EZ MIRROR HANDLE #6 C/S
12. 6950 #4 MIRROR C/S FRONT SURFACE
13. 7423 SHARPENING STONE ARKANSAS
14. 7424 SHARPENING STONE ARKANSAS EDGE
15. 7502N CASSETTE 16 HOLD DOUBLE DECK



Sickle/Jacquette Scalers  
Black Titaniums  
Sickle Scalers/Curettes  
Scaler Curettes  
Universal Curettes  
Gracey Curettes  
Debridement/Furcations  
Sharpening Stones

HANDLE OPTION

\* 6EZ D/E OR S/E

To remove calculus on the surface of gingival area.



To remove calculus on the surface of gingival area.



1201  
H6/7



1223  
N137  
ANTERIOR



1224  
M0/00  
ANTERIOR



1227  
NV-1/H5



1228  
MTJC



1225  
NV-1



1226  
NV-2

## BLACK TITANIUM INSTRUMENTS

Our black titanium instruments are enhanced in strength and durability up to 3 times more than the ordinary tips. It is highly hygienic as it protects the instruments from virus infection. It is easy to clean and requires less sharpening than other instruments.

The titanium provides the surface of the instrument to be a special non-stick design.



1601T  
GR1/2  
BLACK TITANIUM

## BLACK TITANIUMS



1201T  
H6/7  
BLACK TITANIUM



1203T  
H5/33  
BLACK TITANIUM



1208T  
204S  
BLACK TITANIUM



1226T  
NV-2  
BLACK TITANIUM



1322T  
MC13/14S  
BLACK TITANIUM



1415T  
COL4R/4L  
BLACK TITANIUM



1416T  
COL13/14  
BLACK TITANIUM



1422T  
BH5/6  
BLACK TITANIUM



1601T  
GR1/2  
BLACK TITANIUM



1606T  
GR11/12  
BLACK TITANIUM



1607T  
GR13/14  
BLACK TITANIUM



1608T  
GR15/16  
BLACK TITANIUM



1609T  
GR17/18  
BLACK TITANIUM



1610T  
AF1/2  
BLACK TITANIUM



1614T  
AF11/12  
BLACK TITANIUM



1615T  
AF13/14  
BLACK TITANIUM



1620T  
SUB1/2  
BLACK TITANIUM



1621T  
SUB11/12  
BLACK TITANIUM



1622T  
SUB3/4  
BLACK TITANIUM



Designed for the removal of supragingival calculus.






SCALER & CURETTES

To remove heavy calculus above or below the surface of gingival area.  
 We angulate and sharpen the blade of our instruments more, to greatly reduce treatment fatigue.



Curettes are types of scalars with which dentists can deliver more precise treatments. The end of the tips are designed with two cutting edges and rounded toes.



LANGER (UNIVERSAL)



PROPHY GRACEY  
SHANKS ARE A LITTLE SHORTER THAN A NORMAL GRACEY



1505  
YG7/8


The end of tips is designed with two cutting edges and rounded toe. It is used on all surfaces

- BUNT: Bunting
- YG: Younger-Good
- BT: Bates
- R: Rule
- MAL: Mallery
- HU: Hutchinson

For the removal of heavy subgingival calculus including deep scaling, root planning and periodontal debridement. It is highly recommended for dentists to be equipped with a full set of Standard Gracey.

**STANDARD GRACEY**


**RIGID & SHORTER  
BLADE AVAILABLE**



**SUB-0 GRACEY**



**AF-GRACEY (AVAILABLE SHORTER BLADE)**








**SUB-0 GRACEY  
SMALL BLADE**





- \* O' Hehir currettes debridement : Removing residual calculus after ultrasonic scaling for deep pockets.
- \* Furcation currettes : Used for curetting the infected areas of periodontal disease.



## India Stone

The coarse grained India Stone is for the initial grinding of very blunt instruments. The shape of the stone enables the grinding of most of the surfaces and radii of the instruments. Fine grinding should then be carried out with a Hart-Arkansas stone.



7423  
SHARPENING STONE  
INDIA SLIP



7424  
SHARPENING STONE  
ARKANSAS FLAT

