


Operative **E**



COMPOSITE SET UP 19 INS (Recommended)

- | | | | |
|-----|-------|-----|-----------------------------|
| 1. | 1006 | 3EZ | EXPLORERS EXP 5 |
| 2. | 2231 | 4EZ | EXCAVATORS 32L |
| 3. | 2501T | 4EZ | CARVERS IPC, TITANIUM |
| 4. | 2702 | 4EZ | CARVING KNIVES 21 |
| 5. | 2814 | 4EZ | PICH |
| 6. | 2901T | 6EZ | PLASTIC FILLING(THIN) 1, TN |
| 7. | 2904T | 6EZ | PLASTIC FILLING(THIN) 4, TN |
| 8. | 2910T | 6EZ | PLASTIC FILLING TUFTS4, TN |
| 9. | 3103 | 6EZ | GINGIVAL RETRACTORS GRM2 |
| 10. | 3108 | 6EZ | CORD PACKER 113 |
| 11. | 4310 | | AMALGAM CARRIERS M/J |
| 12. | 4706 | | FORCEP PAPER |
| 13. | 6517 | | FORCEP LONDON-COLLEGE |
| 14. | 6911 | 6EZ | MIRROR HANDLE #6 C/S |
| 15. | 6951 | | #5 MIRROR C/S FRONT SURFACE |
| 16. | 7213 | | #207 RUBBER DAM CLAMPS |
| 17. | 7217 | | #27 RUBBER DAM CLAMPS |
| 14. | 7233 | | #210 RUBBER DAM CLAMPS |
| 15. | 7236 | | IVORY LIGHT WEIGHT 17CM |
| 16. | 7502 | | CASSETTES 8-16 HOLD |
| 17. | N | | #1 RUBBER DAM CLAMP |
| 18. | N | | #14 RUBBER DAM CLAMP |
| 19. | N | | #9 RUBBER DAM CLAMP |

Chisels
 Excavators
 Plugger/Condensers
 Carvers
 Carving Knives
 Burnishers
 Gingival Retractors
 Cord Packers
 PKT Waxing
 Amalgam Carriers

HANDLE OPTION


- * 4EZ D/E OR S/E
- * 6EZ D/E OR S/E

- * Hatchets : For cavity preparation and removing hard caries.
- * Hoes : For forming line angles on anterior area.
- * Chisels : For refining the cavity preparation.
- * Margin Trimmer : For producing proper bevel on enamel margins.

HATCHETS



CHISELS




MARGIN TRIMMERS



HOE



For removing debris and decay from cavities.



For removing debris and decay from cavities.



2230
EXC 31L



2231
EXC 32L



2232
EXC 33L




2233
EXC 34L




2234
EXC 35L

For packing and filling materials into cavity preparations.



* T A/V : Available Titanium


For packing and filling materials into cavity preparations.




For carving anatomical features, smoothing surfaces and trimming excess materials.




ANDREWS




CLAPPISON



HOLLENBACK CLEOID




HOLLENBACK




SHOSHAN



THARP



For trimming excess filling materials.





For condensing and smoothing carved amalgam.




■ M: Mini ■ B: Boghosian ■ Titanium Available

For composite placing and shaping.

THIN COMPOSITE



UNIVERSAL




- M: Mini ■ B: Boghosian ■ T/G: Titanium–Gold Available
- B: BEQ(Quetin) ■ T: Tufts ■ WDS: Woodson ■ Titanium–Available

For composite placing and shaping.

ANTERIOR



POSTERIOR



- * Gingival Retractors : For protecting tissue during cavity preparations and placement of composite fillings.
- * Cord Packer : For protecting tissue of gingival area by putting the thread between teeth and gingival area.

GINGIVAL RETRACTORS



CORD PACKERS



For shaping, carving and flowing of molten wax.



3203
3R



3201
1



3202
2



3204
3



3205
4



3206
5



3211
GOLD FOGEL-1

For carrying and dispensing amalgam materials.



4301
SMALL 1.5mm



4302
MEDIUM 2mm



4303
LARGE 2.5mm



4304
JUMBO 3mm



4305 **4306** **4307**

REMOVABLE TIP



4311
AMALGAM WELL SMALL



4308
-0.15mm- SMALL/MEDIUM -2mm-



4310
-2mm- MEDIUM/JUMBO -2.5mm-



4310 EZ
-2mm- EZ MEDIUM/JUMBO -2.5mm-



2404T
PLG-3 TITANIUM



2407T
PLG-3 TITANIUM



2413T
PLGO-1 TITANIUM



2415T
PLGO-3 TITANIUM



2415TR
PLGO-TR TITANIUM



2420T
PLG S/M TITANIUM



2421T
PLG M/L TITANIUM



2422T
KTO TITANIUM



2501T
IPC TITANIUM



2502T
IPCOA TITANIUM



2801T
B-18 TITANIUM



2802T
B26/27S TITANIUM



2803T
B27/29 TITANIUM



2806T
B21 TITANIUM



2807T
B21L TITANIUM



2820T
B28/29 TITANIUM



2821T
B29/34 TITANIUM



2825T
RU-1 TITANIUM



2826T
RU-2 TITANIUM



2827T
RU-3 TITANIUM



2901T
PFIT-1 TITANIUM



2902T
PFIT-2 TITANIUM



2903T
PFIT-3 TITANIUM



2904T
PFIT-4 TITANIUM



2910T
PFI-4 TITANIUM



2920T
WDS-2 TITANIUM



2921T
PFI CM-11 TITANIUM



2922T
PFIM-12 TITANIUM



2923T
PFIT-5 TITANIUM



3001T
PF18A TITANIUM



3002T
PFIT-11 TITANIUM